

# Clear Surgical Tape



Sterex Medical plastic tape is a transparent polyethylene film tape with strong adhesion and a bi-directional tear. It is perforated to promote greater skin breathability and is hypoallergenic.

Since It allows for visualization of the skin and is very easy to tear into thin strips, Sterex transparent tape is ideal for use on I.V. or catheter sites.

Product Code	Unit of Measure	Units per Box	Boxes per Case	NDC #
TRNSP1	1"	12	10	
TRNSP2	2"	12	10	

**Active Ingredient:**