

# Sterile Alcohol Prep Pads



Sterex Sterile Alcohol Prep Pads act as a topical anti-septic. They are ideal for cleansing and disinfecting skin and in the preparation of site for insertion or incision. The waffle style indentations on the 1" x 1 1/4" non-woven pads enhance surface friction for added cleansing and protection against infection. Each 2-ply cloth-like wipe is saturated with 70% Isopropyl alcohol.

**All Sterex Medical Alcohol Preps come sealed in airtight, poly-lined foil pouches to ensure sterility and to minimize dried out pads**

Product Code	Unit of Measure	Units per Box	Boxes per Case	NDC #
ALCPS1	EA	200/BX	10BX/CS	N/A

## ***Active Ingredients:***

### **Isopropyl Alcohol**

Isopropyl Alcohol is widely used by medical professionals as a dermal antiseptic and disinfectant and can kill many different kinds of bacteria, including some antibiotic resistant bacteria and TB bacteria. A 70% mixture of isopropyl alcohol can kill non-spore forming bacteria, virus and fungus in as fast as 15-30 seconds both on surfaces and on human skin.